

Cardioprotective Haemodialysis

# Media Supply

At your fingertips



Cardioprotective Haemodialysis **SPOT**




**FRESENIUS  
MEDICAL CARE**

# Protect your Patient

## Media Supply – the all-in-one distribution system

The Online Purification Cascade (OPC) is a significant part of Fresenius Medical Care's integrated approach to dialysis therapy. The OPC product line covers all processes for the delivery of high-quality fluids required for dialysis treatment.

Media Supply products provide the adequate interface for easy access to all dialysis-relevant media like purified dialysis fluid, concentrates or electricity.



**For dialysis, a hygienic and orderly treatment area is essential. Fresenius Medical Care's integrated Media Supply solution delivers dialysis media to the point of care in a clean and organised manner.**

### Convenient


- A covered drain reduces hygienic risks, noise and unwanted odours
- Quick and easy access to all interfaces
- Easy to install, extend or modify
- Large choice from a broad portfolio of hydraulic and electrical services

### Protective

- Reduced risk of failure as electric and hydraulic circuits are physically separated
- Free-fall drain prevents microbial back- and cross-contamination
- Easy cleaning and disinfection facilitated by the choice of surface material

### Economic

- Flexible investment facilitated by a wide range of configuration options



**Three Media Supply systems allow for customised installation.**

**SPOT on:**

- Clean and organised delivery of dialysis media at the point of care

### **Rail System**

Media Supply Rail (MSR) is designed for simplicity providing a minimalistic dialysis fluid interface.

---



### **Column System**

Media Supply Centre (MSC) combines the provision of dialysis fluids and electrical services in compact, preconfigured individual columns. It balances small space requirements and maximised performance.

---



### **Panel System**

Permeate, concentrate, drain, alarms, equipotential bonding, Ethernet, lights, audio & TV, and even medical gases are no problem for Media Supply Panel (MSP).

It provides all the functionality you could ask for. Developed for maximum convenience in operation and maintenance, it can be wall-mounted and coloured to blend in with any interior design.

---



## Technical Data

Specifications					
	MSR	MSC		MSP	
<b>Dialysis places</b>	1	1	2	1	2
<b>Permeate connectors</b>	1	1	2	1	2
<b>Concentrate connectors</b>	up to 3	up to 3	2x up to 3	up to 3	2x up to 3
<b>Drain</b>	1	1	2	up to 2	2x up to 2
<b>Dimensions</b> (h x w x d)	42 x 50 x 7 cm	87 x 39 x 19 cm	87 x 51 x 19 cm	88 x 51 x 19 cm	88 x 51 x 19 cm
<b>Weight</b>	max. 2 kg	max. 15 kg	max. 15 kg	max. 15 kg	max. 15 kg

If necessary, the exact configuration of electrical and hydraulic services can be adapted individually.

